



ASPARAGINASE REQUISITION FORM
601 Biotech Drive, Suite 301
North Chesterfield, VA 23235
PHONE: (844)347-2643; FAX: (804)977-5041
Email: clientservices@grangergenetics.com

**** SPECIMEN REQUIREMENTS**

Specimen: Serum

Collection: Venous Collection, Red-top tube or gel-barrier tube, separate serum from RBCs within 2 hours, Freeze/Refrigerate serum immediately

Volume: 0.5 – 2 mL

Storage Instructions: Freeze prior to transport. Sample may be refrigerated at 2°C to 8°C, see stability below.

Temperature	Timeframe (activity assay)	Timeframe (antibody assay)
Room Temperature	1 Day	1 Day
Refrigerated	3 Days	3 Days
Frozen	30 Days	30 Days
Freeze/Thaw Cycles	3 Cycles	3 Cycles

*****SHIPPING RECOMENDATION**

Sample Temperature: It's recommended to send the samples frozen but cold is fine if temperature can be maintained during transit.

Shipping: Samples should be sent via Priority Overnight Express delivery. **First AM shipping is not necessary.**

- Shipping on dry ice is preferred. If using cold packs, please freeze sample prior to shipping.

Sample Delivery: Ship samples to be delivered Monday through Friday.

- ***PLEASE NOTE: Granger Genetics is closed and does not receive samples on Saturday or Sunday, please plan your shipments accordingly.***
- ***Any additional closings, such as Holidays, will be posted on our website.***